

IVUS

Delivering more clinically relevant information

IVUS

Delivering more clinically relevant information

Why IVUS?

IVUS

Delivering more clinically relevant information

*Fractional Flow Reserve, BSC Data on file

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

- 1 Abizaid, Andrea et al, J AM Coll Cardiol 1999; 34:707-715
- 2 Mintz, Intracoronary Ultrasound, 2005;122
- 3 Rizik, et al Benefits of Cutting Balloon Before Stenting *JINVAS CARDIOL* 2003;15:624-628
- 4 Moussa, Moses, Colombo et al. Coronary Stenting After Rotational Atherectomy in Calcified and Complex Lesions. *Circulation* 1997; 96: 128-136
- 5 Costa, M, et. al *American Journal of Cardiology* 2008; 12: 1704-1711
- 6 CRUISE Study, 2003
- 7 Chenau et al, 2003
- 8 Broddie B, Cooper C, Jones M, Fitzgerald P, Cummins F, et al. Is adjunctive balloon post-dilatation necessary after coronary stent deployment? Final results from the POSTIT Trial
- 9 Roy et al. Eur Heart J 2008; 29:1851-7
- 10 Costantini et al. TCT 2008 Abstract
- 11 Mintz, Gary, TCT 2008, IVUS Use in the DES Era, Routine or Selective Use
- 12 Hermiller, James, MD, May 2008 Case Study

European Distribution Centre –

The Netherlands T: +31 45 54 67 700 F: +31 45 54 67 800

European Headquarters –

Paris T: +33 1 57 66 80 00 F: +33 1 57 66 84 99

Argentina (Freephone) T: 0800 555 2878 F: +5411 4896 8517

Australia/New Zealand T: +61 2 8063 9100 F: +61 2 9330 1404

Austria T: +43 1 60 810 F: +43 1 60 810 60

Belgium (Freephone) T: 0800 94 494 F: 0800 93 343

Brazil T: +55 11 5502 8500 F: +55 11 5103 2212

Canada T: +1 888 359 9691 F: +1 888 575 7396

Chile T: +562 445 4904 F: +562 445 4915

China – Beijing T: +86 10 8525 1588 F: +86 10 8525 1566

China – Guangzhou T: +86 20 8767 9791 F: +86 20 8767 9789

China – Shanghai T: +86 21 6141 5959 F: +86 21 6141 5900

Colombia T: +57 1 629 5045 F: +57 1 612 4761

Czech Republic T: +420 296 331 901 F: +420 296 331 935

Denmark (Freephone) T: 80 30 80 02 F: 80 30 80 05

Eire (Freephone) T: 1800 882 969 F: 1800 882 968

Finland T: +358 20 762 8882 F: +358 20 762 8883

France T: +33 1 39 30 49 00 F: +33 1 39 30 49 01

Germany (Freephone) T: 0800 723300 F: 08000 723319

Greece T: +30 210 9542 300 F: +30 210 9542 310

Hong Kong T: +852 2960 7100 F: +852 2563 5276

Hungary T: +36 1 456 30 40 F: +36 1 456 30 41

India – Bangalore T: +91 80 2212 4928/9 F: +91 80 2207 5153

India – Chennai T: +91 44 2220 1879 F: +91 44 2220 0716

India – Delhi T: +91 11 4243 2222 F: +91 11 2610 0808

India – Mumbai T: +91 22 4030 9169 F: +91 22 4040 9199

Italy T: +39 010 60 60 1 F: +39 010 60 60 200

Korea T: +82 2 3476 2121 F: +82 2 3476 1776

Mexico T: +52 55 5687 63 90 F: +52 55 5687 62 28

Middle East/Gulf/North Africa T: +961 1 805 410 F: +961 1 805 445

The Netherlands T: +31 30 602 55 44 F: +31 30 602 55 05

Norway (Freephone) T: 800 104 04 F: 800 101 90

Poland T: +48 22 435 14 14 F: +48 22 435 14 10

Portugal T: +351 1 381 25 40 F: +351 21 381 25 58

South Africa T: +27 11 640 1600 F: +27 11 463 6077

South East Asia – Malaysia T: +60 3 2283 3813 F: +60 3 2284 3813

South East Asia – Philippines T: +63 2 687 8994 F: +63 2 687 3047

South East Asia – Singapore T: +65 6418 8888 F: +65 6418 8889

South East Asia – Thailand T: +66 2 6390 100 F: +66 2 6390 400

South East Europe T: +30 210 9542 300 F: +30 210 9542 310

Spain T: +34 901 11 12 15 F: +34 91 319 50 03

Sweden T: +46 42 25 69 00 F: +46 42 25 69 69

Switzerland T: +41 32 626 57 00 F: +41 32 626 57 01

Taiwan T: +886 2 2747 7278 F: +886 2 2747 7270

Turkey T: +90 216 464 36 66 F: +90 216 464 36 67

UK (Freephone) T: 0844 800 4512 F: 0844 800 4513


Uruguay T: +59 82 900 6212 F: +59 82 900 6212

Venezuela T: +58 212 959 6275 F: +58 212 959 5328

Boston
Scientific

Defining tomorrow, today.™

www.bostonscientific-international.com




Reasons to Use IVUS

IVUS as a Diagnostic Tool

Angiography provides a two-dimensional image whereas IVUS offers a three-dimensional detailed image with additional information that can:



- Evaluate lesion characteristics - calcification, fibrous, lipid/necrotic core, mixed plaque, eccentricity, size, length
- Impact treatment choices, including:
 - Pre-dilatation
 - Plaque modification
 - Stent selection, size and length
 - Medical management
 - No treatment



Results from case studies are not predictive of results in other cases.
Results in other cases may vary.

Assessment of Lesion Significance

Comparing IVUS to FFR and based on studies with a clinical outcome, MLA less than 4.0mm² of a proximal epicardial artery is considered to be a flow-limiting stenosis.



Results from case studies are not predictive of results in other cases.
Results in other cases may vary.

IVUS empowers you to treat the patient not the number.

IVUS Gives You More		
Measurement	FFR	IVUS
Lesion Significance	X	X
Stent Expansion		X
Stent Placement		X
Stent Apposition		X
Stent Diameter Sizing		X
Stent Length Sizing		X
Plaque Burden		X
Plaque Tissue Type		X
% Stenosis		X
Lesion Morphology		X

Assessment of Left Main

IVUS can be beneficial during these angiographic challenges in the Left Main¹:


- Aortic cusp opacification or streaming of contrast obscures the ostium
- Short length of vessel leaves little normal vessel for comparison
- Distal left main artery may be concealed by the branching of the left anterior descending and the left circumflex

Guidance for Plaque Modification

IVUS helps assess lesion characteristics that may impact treatment choices like³:

- Pre-dilatation
- Pre-dilatation with cutting balloon
- Rotablator™ Rotational Atherectomy System

Lesions which initially appear as either treatable with PTCA or by stenting may benefit from IVUS to help assess plaque morphology. Asymmetrical stent expansion occurs in up to 50% of lesions if plaque modification was not performed before stent deployment.⁴



Results from case studies are not predictive of results in other cases.
Results in other cases may vary.



Guidance for Stenting

IVUS guidance, during stent procedures, is important due to geographic miss, stent apposition and expansion.

- S.T.L.L.R. Trial shows us that 67% of the time there is some kind of geographic miss when landing a stent.⁵
- Studies involving stent apposition show:
 - 44% lower TLR with complete stent apposition⁶
 - 78% of SATs involved incomplete stent apposition⁷
- The POSTIT Trial revealed that when using only the stent delivery balloon, over 70% of patients did not achieve optimal stent deployment⁸



IVUS is Clinically Relevant in DES

Over 3100 patients studied comparing IVUS-guided DES implantation verses angiography alone showed the use of IVUS reduced the need for repeat procedures



Influence on Stent Thrombosis

In two different studies, using IVUS during DES implantation has the potential to influence treatment strategy and affect DES thrombosis.



Assessments of Complex Cases and Complications

IVUS may be beneficial during complex cases or evaluation of complications.

High Risk Patient Subsets¹¹:

- Renal failure
- Limitations to dual antiplatelet therapy use
- Diabetes mellitus
- Poor left ventricular function

High Risk Lesion Subsets¹¹:

- Bifurcations
- Ostial lesions
- Small vessels
- Long lesions
- Treatment of In-stent restenosis

Complications¹²:

- Dissections
- Hematomas



Results from case studies are not predictive of results in other cases.
Results in other cases may vary.