


Shaping Tool

Motor Housing

Control Unit



TruePath™
CTO Device

TruePath™ CTO Device

UPN	Catalog No.	Description
H74939208181650	39208-18165	TruePath CTO Device
H74939208001601	39208-00160	TruePath Extension Wires (5/box)


TruePath™ Compatibility

Size (OD)	Length	Sheath Compatibility	Min. Guide Catheter ID	Tip Profile
0.018" / 0.46mm	165cm	≥ 0.018" / 0.46mm compatible catheters*	≥ 0.018" / 0.46mm	0.017" / 0.43mm

* With usable lengths of ≤ 135cm



Power Forward.



Rotating diamond-coated distal tip.

TruePath™: a powerful new solution for intraluminal treatment of chronic total occlusions. At 0.018"/0.46 mm in diameter, the guidewire-like profile is nearly half the size of competitive offerings, designed for optimal steering and crossing. Once positioned, the self-rotating tip works effortlessly.