

# Gateway™ Advantage Y-Adapter

Boston  
Scientific



# Gateway™ Advantage Y-Adapter

The Gateway Advantage Y-Adapter is designed to prevent irrigant backflow and secure devices inserted through ureterscope working channels.

- Provides a leak-resistant seal around guidewires, laser fibers, retrieval devices and other devices
- Designed for smooth introduction and removal of urological devices
- Y-Adapter design for additional irrigation connection or contrast injection
- No detectable latex



Rotatable male luer connection designed for easy connection to ureterscope port providing ease efficiency and flexibility to adapt to needs of the case

Blue luer seal provides the flexibility of connecting to irrigant tubing or creating a leak-resistant seal around urological devices up to 3F

Tuohy Borst valve with easy-to-adjust thumb screw closure designed for passage of devices up to 8F



Blue luer seal and side port compatible with standard male luer connection

## ORDERING INFORMATION

Order Number	Product Description	Units
M0067301501	Gateway Advantage Y-Adapter	Box 10

Gateway is an unregistered trademark of Boston Scientific Corporation or its affiliates. All other trademarks are property of their respective owners.

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.** Please refer to package insert provided with the product for complete Indications for Use, Contraindications, Warnings, Precautions, Adverse Events, and Instructions prior to using this product.

**Boston Scientific**

Advancing science for life™

Boston Scientific Corporation  
One Boston Scientific Place  
Natick, MA 01760-1537  
[www.bostonscientific.com/urology](http://www.bostonscientific.com/urology)

**Ordering Information**  
**1.888.272.1001**

© 2013 Boston Scientific Corporation or its affiliates. All rights reserved.

URO-151601-AB XM 5/13