

Mach 1TM

Peripheral Guide Catheter



*The
Complete
Performance
Package*

Boston
Scientific

Mach 1™

Peripheral Guide Catheter

The Complete Performance Package.

Convenient handling, delivery and versatility make Mach 1™ the complete performance package providing access, deliverability and visibility that peripheral vascular procedures require.

Precise Handling

- **Proprietary Polymer Technology** is designed to provide curve retention and back-up support
- **Round Double-Braid Design** with tungsten weave aims to provide kink resistance, torque control, pushability and radiopacity
- **Atraumatic Tip** designed to reduce vessel trauma and to allow thromboaspiration technique

Excellent Delivery*


- **Larger Lumens*** provide dye flow for improved visibility
- **Smooth PTFE Liner** features lubricious surface for improved device passage and deliverability

Particular Versatility


- **Broad Size Selection** ranging from 6-8F provides the versatility to address more clinical applications
- **19 Interventional Curve Styles** designed for excellent access and performance of most peripheral procedures
- **Device Compatibility** makes Mach 1™ our guide catheter to choose for a multitude of treatment options

*Competitive Inner lumen diameter data on file

Ordering Information



6 F / 0.070" (OD / ID)			7 F / 0.081" (OD / ID)			8F / 0.091" (OD / ID)		
Order Number	Tip Shape	Length (cm)	Order Number	Tip Shape	Length (cm)	Order Number	Tip Shape	Length (cm)
M001196000	HS	55	M001197000	HS	55	M001198000	HS	55
M001196660	RDC	55	M001197660	RDC	55	M001198660	RDC	55
M001196960	RDC 1	55	M001197960	RDC 1	55	M001198960	RDC 1	55
M001196040	RE-S	55	M001197040	RE-S	55	M001198040	RE-S	55
M001196080	RE-SS	55	M001197080	RE-SS	55	M001198080	RE-SS	55
M001196120	RE-L	55	M001197120	RE-L	55	M001198120	RE-L	55
M001196160	C1	55	M001197160	C1	55	M001198160	C1	55
M001196200	C2	55	M001197200	C2	55	M001198200	C2	55
M001196230	MP1	55	M001197230	MP1	55	M001198230	MP1	55
M001196260	MP1	90	M001197260	MP1	90	M001198260	MP1	90
M001196300	MP2	55	M001197300	MP2	55	M001198300	MP2	55
M001196330	MP2	90	M001197330	MP2	90	M001198330	MP2	90
M001196450	ST	55	M001197450	ST	55	M001198450	ST	55
M001196480	ST	90	M001197480	ST	90	M001198480	ST	90
M001196530	CROSSOVER 1	55	M001197530	CROSSOVER 1	55	M001198530	CROSSOVER 1	55
M001196560	CROSSOVER 2	55	M001197560	CROSSOVER 2	55	M001198560	CROSSOVER 2	55
M001196620	LIMA	55	M001197620	LIMA	55	M001198620	LIMA	55
M001196700	BATES 1	55	M001197700	BATES 1	55	M001198700	BATES 1	55
M001196800	BATES 2	55	M001197800	BATES 2	55	M001198800	BATES 2	55
M001196850	BATES 3	55	M001197850	BATES 3	55	M001198850	BATES 3	55
M001196900	BATES 4	55	M001197900	BATES 4	55	M001198900	BATES 4	55
M001196380	40°	90	M001197380	40°	90	M001198380	40°	90



Boston Scientific

Freephone from Sweden:
Tel.: 020 65 25 30 • Fax: 020 55 25 35

Freephone from Finland:
Tel.: 09 622 42 85 • Fax: 09 622 42 86

Freephone from Denmark:
Tel.: 80 30 80 02 • Fax: 80 30 80 05

Freephone from Norway:
Tel.: 80 010 404 • Fax: 80 010 190

Freephone from Belgium:
Tel.: 080094 494 • Fax: 080093 343

Freephone from Germany:
Tel.: 0800 072 3302 • Fax: 0800 072 3319

European Headquarters (Paris)
Tel: +33 1 46 49 66 00 • Fax: +33 1 46 49 66 99

France
Tel: +33 1 39 30 49 30 • Fax: +33 1 39 30 49 01

Italy
Tel: +39 010 60 60 236 • Fax: +39 010 60 60 200

Switzerland
Tel: +41 32 626 57 00 • Fax: +41 32 626 57 01

Spain
Tel: +34 93 444 7226 • Fax: +34 93 405 9045

Poland
Tel: +48 22 435 14 14 • Fax: +48 22 435 14 10

Czech Republic
Tel: +420 2 3536 2911 • Fax: +420 2 3536 4334

Hungary
Tel: +36 1 456 30 40 • Fax: +36 1 456 30 41

Austria
Tel: +43 1 726 3005 1021 • Fax: +43 1 726 3005 10 15

UK
Tel: +44 1727 866633 • Fax: +44 1727 830371

The Netherlands
Tel: +31 43 356 8260 • Fax: +31 43 356 8265

Portugal
Tel: +351 1 381 25 40 • Fax: +351 1 381 28 29

www.bostonscientific.com

© 2003 Boston Scientific Corporation.
Printed in Germany.
DEUVAS2065EA

™ All trademarks stated in this document are the property of Boston Scientific Corporation.