

5 Ureteral Drainage

Description	Order No.	Order Qty.
Polaris™ Ultra Ureteral Stent 5Fr x 24cm	M0061921220	
Polaris Ultra Ureteral Stent 5Fr x 26cm	M0061921230	
Polaris Ultra Ureteral Stent 6Fr x 24cm	M0061921320	
Polaris Ultra Ureteral Stent 6Fr x 26cm	M0061921330	
Polaris Ultra Ureteral Stent 7Fr x 24cm	M0061921420	
Polaris Ultra Ureteral Stent 7Fr x 26cm	M0061921430	
Polaris Ultra Ureteral Stent 8Fr x 24cm	M0061921520	
Polaris Ultra Ureteral Stent 8Fr x 26cm	M0061921530	

6 Accessories

Description	Order No.	Order Qty.
SAPSTM (Single Action Pumping System) (Box5)	M0067201001	
SAPS (Single Action Pumping System) Continuous Flow (Box 5)	M0067201011	
Dual Lumen Catheter 10Fr X 50cm	M0064051000	
Urolok™ II Adaptor (Box 10) w/ Ball-Fitting Adaptor	M0067301401	
Encore® Inflation Device (Box 5) 20cc w/ Pressure Gauge	M0067101131	
LeVeen™ Inflation Device (Box 5) 10cc w/ Pressure Gauge	M0067101111	

HOSPITAL ORDER FORM

Physician Name: _____

Physician Signature: _____

Hospital Name: _____

Date: _____

Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/urology

Ordering Information
1.888.272.1001

© 2005 Boston Scientific Corporation or its affiliates. All rights reserved.

DVU1030 Rev. A 7.5M 05/05-05/07

Glidewire Guidewire manufactured by Terumo Corporation and distributed by Boston Scientific Corporation.

Sharplan, SlimLine, SlimLine EZ are manufactured by Lumenis and exclusively distributed by Boston Scientific Corporation.

AccuFlex Fibers are manufactured by InnovaQuartz and exclusively distributed by Boston Scientific Corporation in the U.S.

Swiss LithoClast FlexProbe is manufactured by EMS and distributed by Boston Scientific Corporation.

Sharplan, Lumenis, SlimLine and SlimLine EZ are trademarks of Lumenis Ltd.

Caution: Federal Law (USA) restricts these devices to sale by or on the order of a physician.

Refer to package insert provided with product for complete instructions for use including Indications for Use, Contraindications, Potential Adverse Effects, Warnings and Precautions prior to using this product.

Educational Programs

PRECEPTORS

The preceptorship program offers physicians the opportunity to view surgery with some of the best skilled key opinion leaders in the urology community. Your local Boston Scientific representative can coordinate your schedule with one of our physician preceptors.

VIEWPOINT IN-SERVICES

Boston Scientific sponsors Viewpoint Techniques in Endourology. Viewpoint is a series of programs designed to educate the operating room staff about less invasive techniques in urology. Many hospitals have found that this program results in improved communication and more efficient procedures. Registered Nurses who attend the Viewpoint program will receive continuing education credits.



Boston Scientific

Delivering what's next.™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/urology

Ordering Information
1.888.272.1001

© 2005 Boston Scientific Corporation or its affiliates. All rights reserved.

DVU1030 Rev. A 7.5M 05/05-05/07

Glidewire Guidewire manufactured by Terumo Corporation and distributed by Boston Scientific Corporation.

Sharplan, SlimLine, SlimLine EZ are manufactured by Lumenis and exclusively distributed by Boston Scientific Corporation.

AccuFlex Fibers are manufactured by InnovaQuartz and exclusively distributed by Boston Scientific Corporation in the U.S.

Swiss LithoClast FlexProbe is manufactured by EMS and distributed by Boston Scientific Corporation.


Sharplan, Lumenis, SlimLine and SlimLine EZ are trademarks of Lumenis Ltd.

Caution: Federal Law (USA) restricts these devices to sale by or on the order of a physician.

Refer to package insert provided with product for complete instructions for use including Indications for Use, Contraindications, Potential Adverse Effects, Warnings and Precautions prior to using this product.

Delivering What's Next™ For Flexible Ureteroscopy

Boston Scientific



1 Initial Access




Sensor® PTFE - Nitinol Guidewire

The **Sensor® Guidewire** and **Stiff Shaft Glidewire® Guidewire** both combine kink-resistant shafts and hydrophilic tips to access difficult anatomy and a stiff shaft to provide exceptional handling characteristics. The **low profile Imager™ II Catheters**, Glidewire Guidewires and Sensor Guidewires, offer a wide line of products with a balance of qualities designed to match particular clinical indications.

Description		Order No.
Imager II Catheter with Straight Tip (Box 5)	5Fr	M006 4003021
Imager II Catheter with Bern Tip (Box 5)	5Fr	M006 4003031
Flexima Open End Ureteral Catheter (Box 20)	5Fr X 70cm	M006 4002011
Stiff Shaft Glidewire Guidewire w/Straight Tip (Box 5)	.025" X 150cm	M006 6301071
Stiff Shaft Glidewire Guidewire w/Angled Tip (Box 5)	.025" X 150cm	M006 6301081
Sensor Dual-Flex Guidewire w/Angled Tip (Box 5)	.035" X 150cm	M006 6703011
Sensor Dual-Flex Guidewire w/Angled Tip (Box 5)	.038" X 150cm	M006 6703021
Sensor Dual-Flex Guidewire w/Straight Tip (Box 5)	.035" X 150cm	M006 6703081
Sensor Dual-Flex Guidewire w/Straight Tip (Box 5)	.038" X 150cm	M006 6703121

2 Dilatation/Access Sheaths



Navigator™ Ureteral Access Sheath

Boston Scientific offers a wide line of dilation products including **ureteral access sheaths**, **balloon dilatation catheters** and **ureteral dilators**. Access Sheaths are the latest addition to this product line, designed to provide easier, less traumatic scope manipulation during ureteroscopy.

Each have been designed to provide an advantageous balance of features. In general, it is necessary to dilate the ureter 1 to 2 French sizes larger than the ureteroscope*.

Even scopes with a very small tip diameter may widen significantly along the body of the scope and require dilatation for access to the proximal ureter.

*Bagley, Demetrius H, Ureteroscopy AUA Update Series, Lesson 26, Vol XIII 1994.

Description		Order No.
Navigator Ureteral Access Sheath (Box 5)	11/13Fr X 36cm	M006 2502031
Navigator Ureteral Access Sheath (Box 5)	13/15Fr X 36cm	M006 2502091
UroMax™ Ultra Balloon Catheter Kit	4cm X 15Fr	M006 2251210
UroMax Ultra Balloon Catheter Kit	4cm X 18Fr	M006 2251220
Passport™ Balloon on a Wire Catheter Kit	2cm X 12Fr	M006 2181200
Passport Balloon on a Wire Catheter Kit	4cm X 12Fr	M006 2181300
Ureteral Dilator Set	6, 8, 10, 12, 14, 16Fr Dilators w/.038" Guidewire	M006 2401000

3 Fragmenting the Stone




AccuFlex™ Single-Use Holmium Laser Fiber

Holmium Lasers have brought stone fragmentation to the forefront of endoscopic stone treatment. Boston Scientific offers a variety of Lithotripsy devices designed to enhance scope flexibility, improve visibility and prevent stone migration.

Boston Scientific and Lumenis (formerly Coherent and ESC/Sharplan) have brought to the market 200µ fibers that offer both flexibility and visibility. The blue **SlimLine 200µ Laser Fiber** is designed for enhanced visibility during treatment of stones in the upper tract. The gold-colored **Sharplan® 200µ Laser Fiber** is designed for flexibility and is a valuable tool for stones in hard-to-reach places. **AccuFlex™ Single-Use Laser Fibers** are designed to offer efficient energy transmission and a high degree of flexibility for use with non-Lumenis lasers.

Description		Order No.
Sharplan/Lumenis Brand Laser Fiber (Gold) (Box 3)	200µ	M006 8408501
SlimLine™ Laser Fiber (Blue) (Each)	200µ	M006 8408400
SlimLine EZ™ Single Use Laser Fiber (Green) (Box 3)	200µ	M006 8408921
AccuFlex™ Single Use Laser Fiber (Box 5)	150µ	M006 8408011
AccuFlex Single Use Laser Fiber (Box 5)	200µ	M006 8408021
Swiss LithoClast® Flex Probe (Each)	2.7Fr X 94cm	M006 8402350
Swiss LithoClast Ultra Flex Probe (Each)	2.7Fr X 94cm	M006 8407390

4 Retrieving the Stone



ZeroTip™ Basket

There have been dramatic advancements in endoscopic instrumentation for the removal of ureteral stones. Boston Scientific has developed a full line of stone retrieval devices that enable the urologist to work effectively.

The **ZeroTip Nitinol Retrieval Basket** features a flat distal surface for atraumatic manipulation and nitinol wire construction offers a kink-resistant, flexible wire design. The Boston Scientific sheath technology combines high strength with a low profile design for maximum scope deflection and irrigation.

The **Escape™ Nitinol Stone Retrieval Basket** is the only retrieval technology on the market that combines both a 1.9Fr tipless nitinol basket with holmium laser compatibility to create one simple retrieval solution.

Description		Order No.
Escape Nitinol Stone Retrieval Basket	1.9Fr X 90cm	M006 3902000
Escape Nitinol Stone Retrieval Basket	1.9Fr X 120cm	M006 3902010
ZeroTip Nitinol Stone Retrieval Basket	1.9Fr X 120cm	M006 3901050
ZeroTip Nitinol Stone Retrieval Basket	2.4Fr X 120cm	M006 3901010
Stone Cone™ Nitinol Retrieval Coil	3.0Fr X 115cm w/ 7mm coil	M006 3903100
Stone Cone Nitinol Retrieval Coil	3.0Fr X 115cm w/ 10mm coil	M006 3903200
Graspit® Nitinol Stone Retrieval Forceps	2.6Fr x 90cm	M006 3204000
Graspit Nitinol Stone Retrieval Forceps	2.6Fr X 120cm	M006 3204010

5 Ureteral Drainage




Polaris™ Ultra Ureteral Stent

Polaris™ Ultra Ureteral Stent is a line extension of the Polaris Ureteral Stent. The Polaris Ultra Stent is made of the same dual durometer Percuflex® Material that provides the physician with the ease of removal benefits of a firmer stent and the bonus of a softer material in the bladder. The Polaris Ultra Stent is designed with an updated Nautilus™ bladder coil to help minimize contact on the sensitive trigone area of the bladder.

*The Polaris Ultra Stent is approved for 365 days indwelling.

Description		Order No.
Polaris Ultra Ureteral Stent	5Fr x 24cm	M006 1921220
Polaris Ultra Ureteral Stent	5Fr x 26cm	M006 1921230
Polaris Ultra Ureteral Stent	6Fr x 24cm	M006 1921320
Polaris Ultra Ureteral Stent	6Fr x 26cm	M006 1921330
Polaris Ultra Ureteral Stent	7Fr x 24cm	M006 1921420
Polaris Ultra Ureteral Stent	7Fr x 26cm	M006 1921430
Polaris Ultra Ureteral Stent	8Fr x 24cm	M006 1921520
Polaris Ultra Ureteral Stent	8Fr x 26cm	M006 1921530

6 Procedure Accessories



SAPS™ Single Action Pumping System

The accessories product line comprises ancillary products which assist in flexible ureteroscopy by providing irrigation and balloon inflation, and helping to combine procedural steps.

The **SAPS (Single Action Pumping System)** is designed to control irrigation (and visualization) with a single hand action clearing pump. The **Dual Lumen Catheter** assists in the advancement a second guidewire, or contrast injection in a single step. The **Encore® Inflation Device** assists in controlled single step balloon inflation, while the **UroLok™ II** scope Adapter helps to smooth device introduction through the flexible ureteroscope. These devices may potentially improve or enhance ureteroscopy.

Description		Order No.
SAPS (Single Action Pumping System) (Box 5)		M006 7201001
SAPS (Single Action Pumping System) Continuous Flow (Box 5)		M006 7201011
Dual Lumen Catheter	10Fr X 50cm	M006 4051000
Urolok II Adaptor (Box 10)	w/ Ball-Fitting Adaptor	M006 7301401
Encore Inflation Device (Box 5)	20cc w/ Pressure Gauge	M006 7101131
LeVeon™ Inflation Device (Box 5)	10cc w/ Pressure Gauge	M006 7101111

1 Initial Access

Description		Order No.	Order Qty.
Imager™ II Catheter with Straight Tip (Box 5)	5Fr	M006 4003021	
Imager II Catheter with Bern Tip (Box 5)	5Fr	M006 4003031	
Flexima Open End Ureteral Catheter (Box 20)	5Fr X 70cm	M006 4002011	
Stiff Shaft Glidewire® Guidewire w/Straight Tip (Box 5)	.025" X 150cm	M006 6301071	
Stiff Shaft Glidewire Guidewire w/Angled Tip (Box 5)	.025" X 150cm	M006 6301081	
Sensor Dual-Flex Guidewire w/Angled Tip (Box 5)	.035" X 150cm	M006 6703011	
Sensor Dual-Flex Guidewire w/Angled Tip (Box 5)	.038" X 150cm	M006 6703021	
Sensor Dual-Flex Guidewire w/Straight Tip (Box 5)	.035" X 150cm	M006 6703081	
Sensor Dual-Flex Guidewire w/Straight Tip (Box 5)	.038" X 150cm	M006 6703121	

2 Dilatation/Access Sheaths

Description		Order No.	Order Qty.
Navigator™ Ureteral Access Sheath (Box 5)	11/13Fr X 36cm	M006 2502031	
Navigator Ureteral Access Sheath (Box 5)	13/15Fr X 36cm	M006 2502091	
UroMax™ Ultra Balloon Catheter Kit	4cm X 15Fr	M006 2251210	
UroMax Ultra Balloon Catheter Kit	4cm X 18Fr	M006 2251220	
Passport™ Balloon on a Wire Catheter Kit	2cm X 12Fr	M006 2181200	
Passport Balloon on a Wire Catheter Kit	4cm X 12Fr	M006 2181300	
Ureteral Dilator Set	6Fr - 16Fr Dilators w/.038" Guidewire	M006 2401000	

3 Fragmenting the Stone

Description		Order No.	Order Qty.
Sharplan/Lumenis Brand Laser Fiber (Box 3)	200µ	M006 8408501	
SlimLine™ Laser Fiber (Each)	200µ	M006 8408400	
SlimLine EZ™ Single Use Laser Fiber (Box 3)	200µ	M006 8408921	
AccuFlex™ Single Use Laser Fiber (Box 5)	150µ	M006 8408011	
AccuFlex Single Use Laser Fiber (Box 5)	200µ	M006 8408021	
Swiss LithoClast® Flex Probe (Each)	2.7Fr X 94cm	M006 8402350	
Swiss LithoClast Ultra Flex Probe (Each)	2.7Fr X 94cm	M006 8407390	

4 Retrieving the Stone

Description		Order No.	Order Qty.
Escape™ Nitinol Stone Retrieval Basket	1.9Fr X 90cm	M006 3902000	
Escape Nitinol Stone Retrieval Basket	1.9Fr X 120cm	M006 3902010	
ZeroTip™ Nitinol Stone Retrieval Basket	1.9Fr X 120cm	M006 3901050	
ZeroTip Nitinol Stone Retrieval Basket	2.4Fr X 120cm	M006 3901010	
Stone Cone™ Nitinol Retrieval Coil	3.0Fr X 115cm with 7mm coil	M006 3903100	
Stone Cone Nitinol Retrieval Coil	3.0Fr X 115cm w/ 10mm coil	M006 3903200	
Graspit® Nitinol Stone Retrieval Forceps	2.6Fr x 90cm	M006 3204000	
Graspit Nitinol Stone Retrieval Forceps	2.6Fr X 120cm	M006 3204010	