



# Stone Cone™

Nitinol Urological Retrieval Coil

Boston  
Scientific



# Stone Cone™

Nitinol Urological Retrieval Coil

The Stone Cone Retrieval Coil is designed for use with intracorporeal lithotripters, including the holmium laser.


The sheath is designed to resist stray holmium laser energy.

Designed to automatically release stones when necessary.

The coil is designed to prevent stone migration during lithotripsy procedures.

Unique design promotes sweeping multiple stone fragments in a single pass.

Quick release handle facilitates scope exchanges.



## Stone Cone Nitinol Retrieval Coil Ordering Information

Product Code	Description
M0063903100	Stone Cone Nitinol Urological Retrieval Coil 3.0Fr, 7mm Coil (Each)
M0063903200	Stone Cone Nitinol Urological Retrieval Coil 3.0Fr, 10mm Coil (Each)

**CAUTION:** Federal Law (USA) restricts this device to sale by or on the order of a physician. Refer to package insert provided with the product for complete Instructions for Use, Contraindications for Use, Potential Adverse Effects, Warnings and Precautions prior to using this product. Stone Cone is a trademark of Med Source Technologies.

# Boston Scientific

*Delivering what's next.™*

Boston Scientific Corporation  
One Boston Scientific Place  
Natick, MA 01760-1537  
[www.bostonscientific.com/urology](http://www.bostonscientific.com/urology)

Ordering Information  
1.888.272.1001

© 2006 Boston Scientific Corporation or its affiliates. All rights reserved.

BVU940 5M 12/06-12/08