




Pulmonary


BIOPSY FORCEPS



NEW Radial Jaw™ 4 Pulmonary Biopsy Forceps




Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515181	Radial Jaw 4 Standard Capacity	1.8	100	2.0	Blue	box 5
M00515182	Radial Jaw 4 Standard Capacity	1.8	100	2.0	Blue	box 20
M00515191	Radial Jaw 4 Standard Capacity w/Needle	1.8	100	2.0	Blue	box 5
M00515192	Radial Jaw 4 Standard Capacity w/Needle	1.8	100	2.0	Blue	box 20
M00515201	Radial Jaw 4 Large Capacity	2.4	100	2.8	Blue	box 5
M00515202	Radial Jaw 4 Large Capacity	2.4	100	2.8	Blue	box 20



Radial Jaw™ 3 Pulmonary Biopsy Forceps




Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515231	Radial Jaw 3	1.8	100	2.0	Blue	box 5
M00515241	Radial Jaw 3 w/Needle	1.8	100	2.0	Blue	box 5
M00515211	Radial Jaw 3 Large Capacity	2.2	100	2.8	Blue	box 5



Radial Edge™ Pulmonary Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00515561	Radial Edge	1.8	100	2.0	Blue	box 5
M00515591	Radial Edge w/Needle	1.8	100	2.0	Blue	box 5




Radial Jaw™ Pulmonary Biopsy Forceps



Order Number	Description	Jaw O.D. (mm)	Working Length (cm)	Minimum Working Channel (mm)	Colour Code	Units
M00512661	Radial Jaw	1.8	100	2.0	Blue	box 5
M00512691	Radial Jaw w/Needle	1.8	100	2.0	Blue	box 5

Pulmonary

TRANSBRONCHIAL ASPIRATION NEEDLES



eXcelon™ Transbronchial Aspiration Needle with Syringe



Order Number	Description	Gauge	Needle		Needle Length (mm)	Catheter Length (cm)	Sheath O.D. (mm)	Units
			O.D. (mm)	I.D. (mm)				
M00564101	eXcelon Transbronchial Aspiration Needle w/Syringe	19	1.07	0.69	15	130	1.8	box 5
M00564111	eXcelon Transbronchial Aspiration Needle w/Syringe	20	0.90	0.58	15	130	1.8	box 5
M00564121	eXcelon Transbronchial Aspiration Needle w/Syringe	21	0.81	0.50	15	130	1.8	box 5

CYTOLOGY BRUSHES


Cellebrity™ Endoscopic Cytology Brushes – Bronchoscope



Order Number	Description	Required Working Channel (mm)	Bristle O.D. (mm)	Sheath Length (cm)	Units
M00516001	Cellebrity Cytology Brush	2.0	1.0	140	box 10
M00516011	Cellebrity Cytology Brush	2.0	1.5	140	box 10
M00516071	Cellebrity Cytology Brush	2.0	1.9	100	box 10
M00516151	Cellebrity Cytology Brush	2.0	1.9	150	box 10

Pulmonary

AIRWAY RETRIEVAL BASKETS



ZeroTip™ Airway Retrieval Baskets



Order Number	Description	Sheath O.D. (mm)	Working Length (cm)	Basket Opening (mm)	Sheath Material
M00513200	ZeroTip Airway Retrieval Basket	0.8	120	12	Polyimide/PTFE
M00513210	ZeroTip Airway Retrieval Basket	1.0	120	16	Polyimide/PTFE

Packaged one per box

PULMONARY DILATION BALLOONS



CRE™ Pulmonary Dilation Balloons




Order Number	Description	Balloon Length (cm)	Initial Balloon O.D. (mm)	Intermediate Balloon O.D. (mm)	Maximum Balloon O.D. (mm)	Working Catheter Length (cm)
M00550300	CRE Balloon	5.5	12	13.5	15	75
M00550310	CRE Balloon	5.5	15	16.5	18	75
M00550320	CRE Balloon	5.5	18	19.0	20	75
M00550330	CRE Balloon	3.0	8	9.0	10	75
M00550340	CRE Balloon	3.0	10	11.0	12	75
M00550350	CRE Balloon	3.0	12	13.5	15	75

Packaged one per box

Pulmonary

INFLATION SYSTEM



Alliance™ II Inflation System

Order Number	Description	Units
M00550620	Alliance™ II Inflation Handle	box 1
M00550601	⊗ 60ml Syringe/Gauge Assembly	box 5

**Airway Stents are located in the Stents Section on pages
33–35 and 37–39**

Pulmonary Guidewires are located in the Guidewires Section on page 48